

Fisher and Paykel SleepStyle Fixed CPAP Machine



The Fisher and Paykel Healthcare SleepStyle Fixed is the latest release in the CPAP industry. The SleepStyle empowers you with the freedom to gain excellent CPAP therapy

We all want the freedom that comes with the ability to choose. To choose how we spend our time. How we go about our work. With The SleepStyle, you can. This machine will be set at your prescribed pressure of 4 – 20cm H₂O

The Fisher and Paykel healthcare Sleep Style Fixed CPAP is powered by technology and strikes the balance between comfort and effective treatment.

The beginning of anything new is usually a crucial time. Learning to adjust to therapy is essential. Fisher and Paykel Healthcare have done everything they can to help, from having central apnea detection, to providing a full range of innovative comfort options.

Quiet, integrated design. The power supply and humidifier have been integrated fully, to minimize space taken on the bedside table. There are no large power packs involved as the small power cable plugs straight to the machine

Get set up and started, simply.

The SleepStyle has simplicity woven into its design. Every detail was considered to make it easy for you to use. In just a few steps, you will be on their way to starting therapy

Easy-access chamber- A simple push of the button gives access to the water chamber which is easy to fill and clean. This is also able to be put in your dishwasher to enhance ease of cleaning.

Ramp - The ramp time and starting ramp pressure are preset to 20 minutes and 1/3 of your fixed pressure requirements but will always start at a minimum of 4cm H₂O

User-friendly menu- In collaboration with current users, Fisher and Paykel Healthcare created a simple menu. It has access to simple changes which you can make to ensure the treatment is more comfortable. It also provides a simple update on the screen every time you turn this off. It will give you an update on the amount of hours used, the Apnea – Hypopnea Index and whether your mask is sealing well by telling you if your mask leakage levels are normal.

ThermoSmart™ - Fisher and Paykel Healthcare are the world leaders in humidification systems. Their world-class humidification features a heated breathing tube with AirSpiral™ technology that is designed to minimize condensation. The enhanced algorithm delivers a comfortable humidification experience for patients. The AirSpiral™ technology is the premier tubing system on the market and provides greater protection against condensation in the tubing and optimises humidification delivery. A simple concept, this tubing is double insulated compared to current heated tubing systems.

SensAwake™ - A responsive ramp throughout the night. When wakefulness is detected, treatment pressure is automatically and promptly reduced to a more comfortable level. This eases the return to sleep and helps effective treatment to continue.

Expiratory Relief -An added layer of comfort. The Expiratory relief automatically relieves CPAP pressure during exhalation to make breathing more comfortable.

Large, Responsive Buttons- The layout and function of the buttons work in harmony with the menu to allow quick navigation and ease of use. All you need to do is just push “Start”.

Data - Get information when you need it. The Fisher and Paykel Healthcare SleepStyle Fixed CPAP has a cellular modem built in which ensures automated data transfer while Bluetooth® technology enables you to track your progress instantly using the SleepStyle App. You’re in control. The patient management and support tools complete a seamless experience to help you succeed in embracing therapy. The machine also has a SmartStick to record all the information needed for your physician..

Keep on track with your CPAP therapy. The Fisher and Paykel healthcare SleepStyle™ App and Web program is designed to enhance engagement with therapy. You will have access to meaningful information allowing you to follow your progress and compare it against how you did on previous nights. You can also take action with troubleshooting tips, making the crucial acclimatization period as easy as possible.

Product Information

Specifications

- Dimensions 144 H x 177 W x 183 D mm (5.6 H x 6.9 W x 7.2 D in)
- Weight 1.7 kg
- Performance Pressure Range: 4 to 20 cmH2O
- Noise Level 28.0 ± 1.5 dBA
- Water Chamber Volume 380 ml up to the maximum water level line
- Data Recording The InfoUSB will store up to 5 years of summary efficacy data, 365 days of detailed efficacy data and 140 hours of high-resolution pressure, leak and flow data. Without an InfoUSB, the device’s internal memory is capable of storing up to 1 year of summary efficacy data, 30 days of detailed efficacy data and 20 hours of high-resolution pressure, leak and flow data.
- Electrical Ratings Rated supply voltage Rated current input Rated supply frequency 100-115 V 1.2 A (2.5 A max.) 50-60 Hz 220-240 V 1.10 A (2.3 A max.) 50-60 Hz

Package Includes

- 1 x Carry Bag 1 x SleepStyle Device
- 1 x ThermoSmart Heated Breathing Tube or Standard Breathing Tube
- 1 x Power Cord 1 x Use and Care Guide
- 1 x Water Chamber
- 1 x Chamber Seal
- 1 x Outlet Seal
- 1 x F&P InfoUSB (already in InfoUSB port)
- 1 x Air Filter (already in device)
- 1 x Spare Air Filter
- 1 x Spare Elbow (for use with a standard breathing tube)
- 1 x SleepStyle Quick Reference Guide